



Metallographic Mounting press Ixiang-5



Description

This machine's optimised heating and cooling system makes the setting time extremely short. And all processes from heating, holding, holding, cooling and unloading are done automatically with one click.

Technology Parameter

Sample Diameter(one size)	φ25mm、φ30mm、φ40mm、φ50mm 1inch,1.25inch,1.5inch,2inch
System Pressure	0-20MPa
Cooling Range	40~90 °C
Heating Temp.	90~200 °C
Heating Power	2400W for 25,30mm 3000W for 40,50mm
Power	4500W
Cooling mode	Water /slow Cooling
Hydraulic press or Mechanical press	Hydraulic press
Stroke	140mm
Voltage	220V, 50Hz
Net weight	44kg
Dimensions	530*485*463mm
Programmable	ten programs can be stored